



Assessment of Key Minerals and Harmful Metals

- Whole Blood, Red Blood Cell and Serum Elements
- Urine Toxic and Essential Elements
- Creatinine Clearance
- Hair Elements
- Fecal Metals



SCIENCE + INSIGHT

Toxic and Essential Elements

Elements are the basic building blocks of all chemical compounds, and human exposure to them occurs both from natural and anthropogenic sources. Many elements are considered nutrients and are essential for the proper functioning of the body. These are generally divided between macrominerals such as calcium, magnesium, potassium, sodium and zinc, and trace minerals including selenium, iodine, boron and molybdenum.

Conversely, there are a number of elements that are toxic to the human body, interfere with its functioning and undermine health—such as mercury, lead, cadmium, aluminum, and arsenic. These toxic metals have no known physiological functions. They can be toxic to organ systems and may disrupt the balance of essential nutrients. Toxic metals and essential element status can be assessed in urine, blood, feces and hair.

Doctor's Data has always employed the best-available techniques as a specialist and pioneer in essential and toxic elemental testing. In fact, we were one of the first clinical reference laboratories in the world to employ ICP-MS and high-resolution ICP-MS for elemental analysis.

Deficiencies of essential trace elements or excessive amounts of heavy metals in the human body can cause significant health effects.

Comprehensive Blood Elements

The standard for diagnosis of lead, mercury or other metal toxicity or poisoning, whole blood metals are also used to assess recent or ongoing exposure to potentially toxic elements. In addition, blood element analysis is ideal for guiding supplementation, and should be performed before and during metal detoxification to evaluate essential element status to ensure treatment safety and effectiveness.

In addition to whole blood, serum elements are used to assess the status of key elements and electrolytes that have important functions in the extracellular fluid compartment of blood, giving a more complete evaluation of total blood element levels.

Red blood cell (RBC) elements tests are used to assess the status of essential elements with important intracellular functions, such as magnesium, copper and zinc. Deficiencies or excesses of these essential elements affect numerous metabolic processes. RBC element analysis is also useful for the assessment of ongoing or recent exposure to specific toxic metals, such as arsenic, cadmium, lead, methylmercury and thallium, that accumulate preferentially in erythrocytes.

Doctor's Data measures essential and toxic metals using ICP-MS for whole blood and red blood cells, and a highly sensitive and specific chemistry analyzer for serum elements.

	Comprehensive Blood Elements	Whole Blood Elements	Serum Elements	Red Blood Cell Elements
Calcium	✓	✓	✓	✓
Magnesium	✓	✓	✓	✓
Copper	✓	✓		✓
Zinc	✓	✓		✓
Manganese	✓	✓		✓
Lithium	✓	✓		
Chromium	✓	✓		
Selenium	✓	✓		✓
Strontium	✓	✓		
Molybdenum	✓	✓		✓
Sodium	✓		✓	
Potassium	✓		✓	✓
Phosphorus	✓		✓	✓
Iron	✓		✓	✓
Boron				✓
Vanadium	✓	✓		✓
Arsenic	✓	✓		✓
Barium	✓	✓		
Cadmium	✓	✓		✓
Cesium				✓
Cobalt	✓	✓		
Lead	✓	✓		✓
Mercury	✓	✓		✓
Nickel	✓	✓		
Platinum	✓	✓		
Thallium	✓	✓		✓
Tungsten	✓	✓		
Uranium	✓	✓		

Whole Blood and Serum Elements tests are available separately or as part of the Comprehensive Blood Elements profile. Red Blood Cell Elements is available as a separate test.





LAB #: B000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00016
 SEX: Female
 AGE: 34

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Essential Elements; Serum

ESSENTIAL ELEMENTS								
		RESULT/UNIT	REFERENCE INTERVAL	-2SD	-1SD	MEAN	+1SD	+2SD
Calcium	(Ca)	9.1 mg/dL	8.6 - 10.3					
Magnesium	(Mg)	2.0 mg/dL	1.7 - 2.5					
Sodium	(Na)	138 mEq/L	133 - 145					
Potassium	(K)	3.5 mEq/L	3.5 - 5.0					
Phosphorus	(P)	3.7 mg/dL	2.5 - 5.0					
Iron	(Fe)	115 µg/dL	50 - 200					

INFORMATION

Sodium and Potassium
 Sodium (Na⁺) and potassium (K⁺) are electrolytes that affect most metabolic functions. They serve to maintain osmotic pressure and hydration of various body fluid compartments, body pH and regulation of heart and muscle functions. Electrolytes are also involved in oxidation-reduction reactions and participate in essential enzymatic reactions. Electrolytes can be affected by state of hydration. Hemolysis can result in falsely elevated K⁺.

Magnesium
 Magnesium (Mg) is a major intracellular cation that is involved in over three hundred enzymatic reactions in the body. Little is known about the factors affecting serum Mg, but the poor diet/malabsorption, diabetes, hyperthyroidism, alcoholism and diuresis. Increased serum Mg lev Addison's disease.

Calcium
 Although 99% of calcium exists in bones and tee nerve impulses, muscle contraction, coagulation, regulated by parathyroid hormone, and serum C muscle tetany while high Ca levels result in lower Marked variations in serum Ca may result from par kidney disease, and other abnormalities.

Inorganic Phosphorus
 Measurements of serum inorganic phosphorus (p parathyroid gland and kidney diseases, and vit parathyroid hormone, and PO₄ levels are inversely muscle weakness, while elevated PO₄ may be ass

Iron
 Measurements of non-heme, serum iron (Fe) are toxicity and acute or chronic hemochromatosis. T ferritin.

Comments:
 Date Collected: 5/16/2014
 Date Received: 5/17/201
 Date Completed: 5/19/2014

©DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Ave
 0001548



LAB #: B000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00000
 SEX: Female
 DOB: AGE: 68

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Toxic & Essential Elements; Whole Blood

ESSENTIAL AND OTHER ELEMENTS								
		RESULT / UNIT	REFERENCE INTERVAL	PERCENTILE				
				2.5 th	16 th	50 th	84 th	97.5 th
Calcium	(Ca)	5.7 mg/dL	4.8 - 7.1					
Magnesium	(Mg)	3.8 mg/dL	3 - 4.2					
Copper	(Cu)	100 µg/dL	65 - 130					
Zinc	(Zn)	730 µg/dL	480 - 780					
Manganese	(Mn)	7 µg/L	4 - 22					
Chromium	(Cr)	0.21 µg/L	0.2 - 0.8					
Lithium	(Li)	0.45 µg/L	0.4 - 2.0					
Selenium	(Se)	150 µg/L	140 - 350					
Strontium	(Sr)	14 µg/L	10 - 45					
Molybdenum	(Mo)	1.2 µg/L	0.3 - 2.5					
Vanadium	(V)	0.045 µg/L	0.04 - 0.3					

TOXIC METALS				
		RESULT / UNIT	REFERENCE INTERVAL	PERCENTILE
				95 th 99 th
Arsenic	(As)	0.9 µg/L	< 9.0	
Barium	(Ba)	0.8 µg/L	< 4.0	
Cadmium	(Cd)	0.3 µg/L	< 1.0	
Cobalt	(Co)	0.1 µg/L	< 0.8	
Lead	(Pb)	7.5 µg/dL	< 3.0	
Mercury	(Hg)	1.7 µg/L	< 4.5	
Nickel	(Ni)	< 1.5 µg/L	< 3.0	
Platinum	(Pt)	< 0.05 µg/L	< 0.10	
Thallium	(Tl)	< 0.05 µg/L	< 0.50	
Tungsten	(W)	< 0.03 µg/L	< 0.10	
Uranium	(U)	< 0.02 µg/L	< 0.10	

Comments:
 Date Collected: 02/02/2017 Time Collected: 1:00 PM Methodology: ICP-MS
 Date Received: 02/07/2017 Fasting: Random
 Date Reported: 02/09/2017
 Blood lead levels in the range of 5-9 µg/dL have been associated with adverse health effects in children aged 6 years and younger.

©DOCTOR'S DATA, INC. • 3755 Illinois Avenue, St. Charles, IL 60174-2420 USA • LAB DIR: Erio Roth, MD • CLIA ID: 14D0646470
 0002020

Results are presented in a clear, easy-to-understand report which details target ranges and graphically illustrates areas of concern. Result-specific commentary is also provided.



LAB #: B000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00000
 SEX: Female
 DOB:

AGE: 61

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Toxic & Essential Elements; Packed Red Blood Cells

		ESSENTIAL AND OTHER ELEMENTS		PERCENTILE				
		RESULT / UNIT	REFERENCE INTERVAL	2.5 th	16 th	50 th	84 th	97.5 th
Calcium	(Ca)	11 µg/g	8-26					
Magnesium	(Mg)	40 µg/g	39-59					
Potassium	(K)	84 mEq/L	72-90					
Phosphorus	(P)	628 µg/g	490-670					
Copper	(Cu)	0.678 µg/g	0.52-0.8					
Zinc	(Zn)	8.0 µg/g	7.8-13.8					
Iron	(Fe)	923 µg/g	780-1000					
Manganese	(Mn)	0.013 µg/g	0.009-0.033					
Selenium	(Se)	0.18 µg/g	0.16-0.49					
Boron	(B)	0.057 µg/g	0.01-0.11					
Molybdenum	(Mo)	0.0003 µg/g	0.0002-0.001					

		TOXIC METALS		PERCENTILE	
		RESULT / UNIT	REFERENCE INTERVAL	95 th	99 th
Arsenic	(As)	0.0086 µg/g	< 0.008		
Cadmium	(Cd)	0.0008 µg/g	< 0.002		
Cesium	(Cs)	0.008 µg/g	< 0.01		
Chromium	(Cr)	0.0008 µg/g	< 0.0005		
Lead	(Pb)	0.032 µg/g	< 0.05		
Mercury	(Hg)	0.017 µg/g	< 0.01		
Thallium	(Tl)	0.00007 µg/g	< 0.00005		

SPECIMEN DATA

Comments:

Date Collected: 02/01/2017
 Date Received: 02/07/2017
 Date Reported: 02/07/2017

Methodology: ICP-MS

0002019

©DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Avenue, St. Charles, IL 60174-2420 • LAB DIR: Erlo Roth, MD • CLIA ID NO: 14D0646470

For more information about all available tests, clinical information and sample reports, visit doctorsdata.com.

Urine Toxic and Essential Elements

Urine Elements are traditionally used to evaluate exposure to potentially toxic elements and wasting of nutrient elements. Toxic metals do not have any useful physiological function. Instead, they adversely affect virtually every organ system and disrupt the homeostasis of nutrient elements.

Additionally, the comparison of urine element concentrations before and after administration of a chelator can be used to estimate net retention of potentially toxic elements. Subsequent urine element analyses, also following the administration of a chelator, are useful for monitoring the efficacy of metal detoxification therapy. Results are expressed per 24 hours or creatinine corrected to account for urine dilution effects.



LAB #: U000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00001
 SEX: Female
 AGE: 61

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Toxic Metals; Urine

TOXIC METALS				
	RESULT µg/g creat	REFERENCE INTERVAL	WITHIN REFERENCE	OUTSIDE REFERENCE
Aluminum (Al)	210	< 35		
Antimony (Sb)	0.5	< 0.4		
Arsenic (As)	40	< 117		
Barium (Ba)	11	< 7		
Beryllium (Be)	< dl	< 1		
Bismuth (Bi)	0.2	< 15		
Cadmium (Cd)	2.2	< 1		
Cesium (Cs)	8.9	< 10		
Gadolinium (Gd)	0.4	< 0.4		
Lead (Pb)	31	< 2		
Mercury (Hg)	15	< 4		
Nickel (Ni)	22	< 12		
Palladium (Pd)	< dl	< 0.3		
Platinum (Pt)	< dl	< 1		
Tellurium (Te)	< dl	< 0.8		
Thallium (Tl)	0.4	< 0.5		
Thorium (Th)	< dl	< 0.03		
Tin (Sn)	1.9	< 10		
Tungsten (W)	1.2	< 0.4		
Uranium (U)	0.2	< 0.04		

URINE CREATININE							
	RESULT mg/dL	REFERENCE INTERVAL	-2SD	-1SD	MEAN	+1SD	+2SD
Creatinine	26.7	35 - 225					

SPECIMEN DATA		
Comments:		
Date Collected: 5/16/2014	pH upon receipt: Acceptable	Collection Period: timed: 6 hours
Date Received: 5/17/2014	<dl: less than detection limit	Volume:
Date Completed: 5/19/2014	Provoking Agent: DMPS CAEDTA	Provocation: POST PROVOCATIVE
Method: ICP-MS	Creatinine by Jaffe Method	
Results are creatinine corrected to account for urine dilution variations. Reference intervals and corresponding graphs are representative of a healthy population under non-provoked conditions. Chelation (provocation) agents can increase urinary excretion of metals/elements.		
V13		



LAB #: U000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00001
 SEX: Female
 AGE: 61

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Essential Elements; Urine

ESSENTIAL AND OTHER ELEMENTS						
	RESULT/UNIT	REFERENCE	PERCENTILE			
	per creatinine	INTERVAL	2.5 th	16 th	50 th	84 th 97.5 th
Sodium (Na)	330 mEq/g	43.5 - 226				
Potassium (K)	79 mEq/g	22 - 82				
Phosphorus (P)	530 µg/mg	250 - 1300				
Calcium (Ca)	1040 µg/mg	35 - 350				
Magnesium (Mg)	480 µg/mg	25 - 230				
Zinc (Zn)	34 µg/mg	0.1 - 2				
Copper (Cu)	0.6 µg/mg	0.01 - 0.09				
Sulfur (S)	1490 µg/mg	308 - 1650				
Manganese (Mn)	0.099 µg/mg	0.0005 - 0.01				
Molybdenum (Mo)	0.12 µg/mg	0.016 - 0.18				
Boron (B)	1.3 µg/mg	0.8 - 6.8				
Chromium (Cr)	0.003 µg/mg	0.0005 - 0.01				
Lithium (Li)	0.023 µg/mg	0.01 - 0.2				
Selenium (Se)	0.18 µg/mg	0.034 - 0.28				
Strontium (Sr)	0.41 µg/mg	0.06 - 0.54				
Vanadium (V)	0.002 µg/mg	0.0002 - 0.004				
Cobalt (Co)	1.9 µg/mg	< 0.008			68 th	95 th
Iron (Fe)	3 µg/mg	< 2				

URINE CREATININE					
	RESULT	REFERENCE	PERCENTILE		
	mg/dL	INTERVAL	-2SD	-1SD	MEAN +1SD +2SD
Creatinine	26.7	35 - 225			

SPECIMEN DATA

Comments:

Date Collected: 5/16/2014 pH Upon Receipt: Acceptable Collection Period: timed: 6 hours
 Date Received: 5/17/2014 <d: less than detection limit Volume:
 Date Completed: 5/19/2014 Provoking Agent: DMPS CAEDTA Provocation: POST PROVOCATIVE
 Method: ISE; Na, K Spectrophotometry; P ICP-MS; B, Ca, Cr, Co, Cu, Fe, Mg, Mn, Mo, Se, Sr, S, V, Zn Creatinine by Jaffe method

Results are creatinine corrected to account for urine dilution variations. Reference intervals and corresponding graphs are representative of a healthy population under non-provoked conditions. Chelation (provocation) agents can increase urinary excretion of metals/elements. V13

DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Avenue, St. Charles, IL 60174-2420 • CLIA ID NO: 14D0646470 • MEDICARE PROVIDER NO: 148453
 0001542

Results are presented in a clear, easy-to-understand report which graphically illustrates target ranges and areas of concern. Result-specific commentary is provided.

For more information about all available tests, clinical information and sample reports, visit doctorsdata.com.

Creatinine Clearance

The Creatinine Clearance test is the most widely used test for estimating glomerular filtration rate (GFR) and renal function. GFR assessment is highly recommended for weighing the advisability of prescribing a variety of drugs, including chelating agents.

The Creatinine Clearance test analyses creatinine in a timed urine collection and a single serum specimen collected during the same period.

Results are presented in a clear, easy-to-understand report complete with result-specific commentary.



LAB #: U000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT -S-00716
 SEX: Female
 AGE: 42

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Creatinine Clearance

RESULTS							
	RESULT / UNIT	REFERENCE INTERVAL	-2SD	-1SD	MEAN	+1SD	+2SD
Creatinine Clearance	118 mL/min	75 - 120					
Urine Creatinine	1300 mg/time	600 - 1900					
Serum Creatinine	0.77 mg/dL	0.6 - 1.3					

INFORMATION

Creatinine Clearance is the most widely used test for estimating glomerular filtration rate (GFR). Creatinine, a breakdown product of muscle creatine, is present in relatively stable levels in serum. It is filtered by the glomeruli and not reabsorbed by the tubules. Changes in renal function are reflected in levels of serum urea and creatinine.

It is not uncommon for elderly patients, and those with heavy metal toxicity to have mild to moderate impairment of renal function. Renal disease is asymptomatic in most cases until late in its clinical course. Safe chelation therapy is highly dependent upon the adequacy of renal function. Excessive mobilization of toxic metals to poorly functioning kidneys may result in renal complications. It is advised that creatinine clearance be monitored prior to and throughout chelation therapy.

Interpretive guidelines:

- 100 mL/min or higher usually indicates normal renal function.
- 50 mL/min or below is indicative of impaired kidney function.
- 30 mL/min or below is indicative of symptomatic renal failure.

Exercise may cause increased clearance. Inaccurate results may be caused by failure to accurately follow the specimen collection instructions.

The calculation for corrected creatinine clearance in mL/min: =

Urine volume per minute x urine creatinine ÷ Serum creatinine x 1.73/body surface area

References:

- Kaplan, Lawrence A., Clinical Chemistry, 3rd Edition. Mosby, St. Louis, 1996
- Jacobs, D.S., Laboratory Test Handbook. 2nd Edition. Lexi-Comp Inc. 1990

SPECIMEN DATA

Comments:

Date Collected: 5/16/2014
 Date Received: 5/17/2014
 Date Completed: 5/19/2014

Height: 0 in
 Weight: 0 lbs
 Body Surface Area: 1.73

Collection Period: 24 hours
 Volume: 900 ml
 Methodology: Automated Jaffe

©DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Avenue, St. Charles, IL 60174-2420 • CLIA ID NO: 14D0646470 • MEDICARE PROVIDER NO: 148453

0001548

For more information about all available tests, clinical information and sample reports, visit doctorsdata.com.

Hair Elements

Hair Elements analysis provides information regarding recent and ongoing exposure to potentially toxic metals, especially methylmercury and arsenic, and time-averaged status of specific nutrient elements. This noninvasive screening test requires only .25 grams of hair. Doctor's Data offers a Hair Elements profile containing essential and toxic elements and a Hair Toxic Element Exposure profile containing an expanded lineup of toxic metals.

A specialist and pioneer in essential and toxic elemental testing since 1972, Doctor's Data has been validated as a supplier of trace element results for the certification of a hair reference material to the European Commission Joint Research Centre.

This noninvasive screening test requires just .25 grams of hair. Choose between the original Hair Elements profile and the expanded Hair Toxic Element Exposure profile. Results are clearly presented along with result-specific commentary.



LAB #: H000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00001
 SEX: Female
 AGE: 51

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Toxic Element Exposure Profile; Hair

TOXIC METALS			PERCENTILE	
	RESULT µg/g	REFERENCE INTERVAL	68 th	95 th
Arsenic (As)	0.021	< 0.14		
Lead (Pb)	0.38	< 3.0		
Mercury (Hg)	0.21	< 3.0		
Cadmium (Cd)	0.032	< 0.20		
Chromium (Cr)	0.52	< 2.0		
Beryllium (Be)	< 0.01	< 0.020		
Cobalt (Co)	0.010	< 0.020		
Nickel (Ni)	0.54	< 2.0		
Zinc (Zn)	170	< 200		
Copper (Cu)	160	< 200		
Thorium (Th)	< 0.001	< 0.002		
Thallium (Tl)	< 0.001	< 0.002		
Barium (Ba)	1.3	< 2.0		
Cesium (Cs)	< 0.002	< 0.005		
Manganese (Mn)	0.19	< 0.50		
Selenium (Se)	0.70	< 1.0		
Bismuth (Bi)	0.018	< 0.020		
Vanadium (V)	0.049	< 0.10		
Silver (Ag)	0.86	< 1.0		
Antimony (Sb)	< 0.01	< 0.020		
Palladium (Pd)	0.011	< 0.020		
Aluminum (Al)	24	< 50		
Platinum (Pt)	< 0.003	< 0.005		
Tungsten (W)	< 0.001	< 0.002		
Tin (Sn)	0.38	< 1.0		
Uranium (U)	0.26	< 0.50		
Gold (Au)	0.082	< 0.10		
Tellurium (Te)	< 0.05	< 0.10		
Germanium (Ge)	0.029	< 0.050		
Titanium (Ti)	0.70	< 1.0		
Gadolinium (Gd)	< 0.001	< 0.002		

SPECIMEN DATA

Comments:
 Date Collected: 5/16/2014 Method: ICP-MS
 Date Received: 5/17/2014 <d: less than detection limit
 Date Completed: 5/19/2014 µg/g = ppm
 Metals are listed in descending priority order based upon data from the Agency for Toxic Substances and Disease
 which considers not only the relative toxicity per gram metal, but also the frequency for occurrence of exposure.

©DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Avenue, St. Charles, IL 60174-2420 • CLIA



LAB #: H000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00001
 SEX: Male
 AGE: 9

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Toxic & Essential Elements; Hair

TOXIC METALS			PERCENTILE	
	RESULT µg/g	REFERENCE INTERVAL	68 th	95 th
Aluminum (Al)	9.0	< 8.0		
Antimony (Sb)	0.088	< 0.066		
Arsenic (As)	0.14	< 0.080		
Barium (Ba)	0.30	< 0.75		
Beryllium (Be)	< 0.01	< 0.020		
Bismuth (Bi)	0.13	< 2.0		
Cadmium (Cd)	0.025	< 0.070		
Lead (Pb)	0.92	< 1.0		
Mercury (Hg)	1.1	< 0.40		
Platinum (Pt)	< 0.003	< 0.005		
Thallium (Tl)	< 0.001	< 0.002		
Thorium (Th)	< 0.001	< 0.002		
Uranium (U)	0.010	< 0.060		
Nickel (Ni)	0.13	< 0.20		
Silver (Ag)	0.14	< 0.14		
Tin (Sn)	0.32	< 0.30		
Titanium (Ti)	0.51	< 0.70		

ESSENTIAL AND OTHER ELEMENTS			PERCENTILE				
	RESULT µg/g	REFERENCE INTERVAL	2.5 th	16 th	50 th	84 th	97.5 th
Calcium (Ca)	157	160- 500					
Magnesium (Mg)	11	12- 50					
Sodium (Na)	100	20- 200					
Potassium (K)	100	12- 140					
Copper (Cu)	11	11- 32					
Zinc (Zn)	350	110- 190					
Manganese (Mn)	0.28	0.08- 0.50					
Chromium (Cr)	0.60	0.40- 0.70					
Vanadium (V)	0.079	0.025- 0.10					
Molybdenum (Mo)	0.14	0.040- 0.090					
Boron (B)	3.6	0.50- 3.5					
Iodine (I)	0.48	0.25- 1.3					
Lithium (Li)	0.010	0.007- 0.020					
Phosphorus (P)	146	150- 220					
Selenium (Se)	0.84	0.70- 1.1					
Strontium (Sr)	0.21	0.21- 2.1					
Sulfur (S)	50900	44000- 51000					
Cobalt (Co)	0.009	0.004- 0.020					
Iron (Fe)	10	7.0- 16					
Germanium (Ge)	0.028	0.030- 0.040					
Rubidium (Rb)	0.086	0.008- 0.080					
Zirconium (Zr)	0.42	0.060- 0.70					

SPECIMEN DATA		RATIOS	
COMMENTS:			
Date Collected: 5/16/2014	Sample Size: 0.198 g	Ca/Mg	14.3
Date Received: 5/17/2014	Sample Type: Head	Ca/P	1.08
Date Completed: 5/19/2014	Hair Color: Brown	Na/K	1
Methodology: ICP/MS	Treatment: Shampoo	Zn/Cu	31.8
		Zn/Cd	> 999

©DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Avenue, St. Charles, IL 60174-2420 • CLIA ID NO: 14D0646470 • MEDICARE PROVIDER NO: 148453

Fecal Metals

Fecal elemental analysis provides a direct indication of dietary exposure to toxic metals and indirect information about the potential for toxic metal burden. Chronic, low-level assimilation of toxic metals can result in accumulation in the body. For many toxic metals, fecal (biliary) excretion is the primary natural route of elimination from the body. Elements are measured by ICP-MS and expressed on a dry weight basis to eliminate variability related to water content of the specimen.

Specimen collection is convenient for the patient and requires only a single-step procedure. Results are presented in a clear, easy-to-understand report and includes result-specific commentary.



LAB #: F000000-0000-0
 PATIENT: Sample Patient
 ID: PATIENT-S-00003
 SEX: Male
 AGE: 7

CLIENT #: 12345
 DOCTOR:
 Doctor's Data, Inc.
 3755 Illinois Ave.
 St. Charles, IL 60174

Toxic Metals; Feces

TOXIC METALS				PERCENTILE	
	RESULT mg/kg Dry Wt	REFERENCE INTERVAL		68 th	95 th
Mercury (Hg)	0.031	<.05 w/o amalgams*			
Mercury (Hg)	0.031	<0.5 with amalgams*			
Antimony (Sb)	0.100	< 0.080			
Arsenic (As)	0.20	< 0.30			
Beryllium (Be)	< dl	< 0.009			
Bismuth (Bi)	229.8	< 0.050			
Cadmium (Cd)	0.41	< 0.50			
Copper (Cu)	63	< 60			
Lead (Pb)	0.27	< 0.50			
Nickel (Ni)	11.8	< 8.0			
Platinum (Pt)	< dl	< 0.003			
Thallium (Tl)	0.019	< 0.020			
Tungsten (W)	0.054	< 0.090			
Uranium (U)	0.085	< 0.120			

WATER CONTENT							
	RESULT % H ₂ O	REFERENCE INTERVAL	-2SD	-1SD	MEAN 72.5%	+1SD	+2SD
% Water Content	67.6	60 - 85%					

INFORMATION

Analysis of elements in feces provides a comprehensive evaluation of environmental exposure, accumulation and endogenous detoxification of potentially toxic metals. For several toxic elements such as mercury, cadmium, lead, antimony and uranium, biliary excretion of metals into feces is the primary natural route of elimination from the body. Studies performed at DDI demonstrate that the fecal mercury content and number of amalgam surfaces are highly correlated, as is the case for post-DMPS urine mercury levels and amalgam surface area.

Results are reported as mg/kg dry weight of feces to eliminate the influence of variability in water content of fecal specimens. The reference values that appear in this report have been derived from both published data and in-house studies at DDI. *Due to exposure to mercury in the oral cavity, people with dental amalgams typically have a considerably higher level of mercury in the feces than individuals without dental amalgams; therefore, two reference ranges have been established for mercury.

To provide guidance in interpretation of results, patient values are plotted graphically with respect to percentile distribution of the population base. Since this test reflects both biliary excretion and exposure (metals to which the patient is exposed may not be absorbed), it may not correlate with overt clinical effects. Further testing can assist in determining whether the metals are from endogenous (biliary excretion) or exogenous (oral exposure) sources.

1. Bjorkman, L, Sandborgh-Englund, G, and Ekstrand, J. Mercury in Saliva and Feces after Removal of Amalgam Fillings. *Toxicology & Applied Pharmacology* 144: 156-162 (1997)
2. Zalups, R. Progressive Losses of Renal Mass and the Renal and Hepatic Disposition of Administered Inorganic Mercury. *Toxicology & Applied Pharmacology* 130: 121-131 (1995)
3. Adamiason, E, Piscator, M., and Nogawa, K. Pulmonary and Gastrointestinal Exposure to Cadmium Oxide Dust in a Battery Factory. *Environmental Health Perspectives*, 26: 219-222 (1979)
4. Smith, J., et al., The Kinetics of Intravenously Administered Methyl Mercury in Man. *Toxicology & Applied Pharmacology* 128:251-256 (1994)
5. Bass, D., et al., "Measurement of Mercury in Feces", Poster presentation 1999 AACC

SPECIMEN DATA


Comments:

Date Collected: 5/16/2014 Provocation: Dental Amalgams: **not indicated**
 Date Received: 5/17/2014 Detoxification Agent: Quantity:
 Date Completed: 5/19/2014 Dosage: Methodology: **ICP-MS**

©DOCTOR'S DATA, INC. • ADDRESS: 3755 Illinois Avenue, St. Charles, IL 60174-2420 • CLIA ID NO: 14D0646470 • MEDICARE PROVIDER NO: 148453

0001547

For more information about all available tests, clinical information and sample reports, visit doctorsdata.com.





SCIENCE+INSIGHT

3755 Illinois Avenue
St. Charles, IL 60174-2420

800.323.2784 (US AND CANADA)
0871.218.0052 (UK)
+1.630.377.8139 (GLOBAL)

doctorsdata.com